

focus Magazine Article Quiz 4th Quarter 2018 The Micron Differences By Beatrice Robichaud NBC Approval #35930 CE Broker # 20-635167

1. Today, if you open any dental magazine, you will read about CAD/CAM products.

- a. True
- b. False

2. CAD/CAM products can be scanners, software applications, machines, and even outsourcing services.

- a. True
- b. False

3. One thing that ads and articles have in common when talking about CAD/CAM is they don't talk about units of measurement.

- a. True
- b. False

4. When talking about precision, they will use a unit of measurement that you have heard of: the micron (fflm).

- a. True
- b. False

5. According to a study that Dr. Brenamark did in 1983, a tolerance of +/- 15fflm can be considered a passive fit.

a. True b. False

6. Even when you have a very precise scanner, there are many key points to keep in mind that will affect the precision of your scanner.

a. True

b. False

7. Manufacturers came up with the idea of creating scan bodies because none of the scanners available in the industry can scan the actual implant.

a. True

b. False

8. Your scan bodies need to fit on the actual implant analogs.

- a. True
- b. False

9. In today's industry, 3 software applications are most widely used. There are other options. In this domain, precision **will be affected** by the software you use.

- a. True
- b. False

10. Precision is not just a matter of the numbers. It is a complex concept affected by the equipment you use and the way you use it.

> a. True b. False

Passing quiz grades are worth ½ point documented scientific credit towards your CDT or RG renewal and towards the State of Florida required 18 hours of continuing education. Look for more quizzes in future issues of *focus* magazine!

You must submit your completed quiz to FDLA by fax at (850) 222-3019 or by mail to the address below. For CDT credit, FDLA will report your hours directly to NBC. For State of Florida credit, FDLA will upload the information to CE Broker. Please allow 2 – 4 weeks to process your quiz.

Name:		_ CDT # if applicable:	DOH Lic #:	
Laboratory Name:				
Address/City/State/Zip:				
Phone:	Email: _			

FDLA, 325 John Knox Rd, Ste L103, Tallahassee, Florida 32303 + (850) 224-0711 phone + (850) 222-3019 FAX

This quiz is provided to test the technician's comprehension of the article's content, and does not necessarily serve as an endorsement of the content by FDLA or NBC.